


SPECTROstar® Nano

Ultra-fast Absorbance Spectrometer for Microplates and Cuvettes



**BMG LABTECH**

The Microplate Reader Company



Innovative Absorbance Microplate Reader with Cuvette Port

Ultra-fast determination of
DNA, RNA, ELISA, and more.

The SPECTROstar^{Nano} represents the next evolution in absorbance microplate reader technology.

Equipped with BMG LABTECH's proprietary ultra-fast UV/Vis spectrometer, the SPECTROstar^{Nano} can measure full spectrum absorbance (220 - 1000 nm) in less than 1 second per well in microplates, cuvettes, and low-volume samples.

Unique Technology and New Features

Backed by German engineering and technology, the SPECTROstar^{Nano} is a versatile, automated absorbance reader for microplates, cuvettes, and low-volume samples. New and unique features include:

- Ultra-fast UV/Vis spectrometer
- Measures all microplate formats up to 1536-wells
- Standard cuvette port for individual samples
- Low-volume samples with BMG LABTECH's LVis Plate (2 µL)
- Multi-mode shaking: linear, orbital, and double orbital
- Incubation up to 45°C
- Gas vent
- Robotic compliant software interface
- Powerful MARS Data Analysis Software


Other features include: automated path length correction, well scanning, kinetic readings, assay-specific buttons, and predefined calculation templates for data analysis.

Perfect for All Absorbance Assays

The SPECTROstar^{Nano} is the ideal instrument to read all absorbance assays in a microplate or cuvette. Its rapid full spectrum analysis allows for absorbance assays never before possible on a microplate reader. The most common absorbance assays such as ELISAs, DNA, RNA, protein (Bradford, BCA, Lowry), and beta-galactosidase have predefined quick run assay-specific protocols and buttons. More complex or user specific assays are programmable.

Assay-Specific Buttons

The Control and MARS Data Analysis Software have predefined assay specific protocols. This user-friendly feature enables researchers to carry out experiments with the SPECTROstar^{Nano} quickly and easily. Measurement and analysis are performed with a single mouse click, and the user can assign a new test run or data analysis to a specific assay button.




Predefined assay-specific buttons for one-click measurements

UV/Vis Spectrometer

The SPECTROstar^{Nano} is equipped with BMG LABTECH's proprietary UV/Vis spectrometer. This spectrometer instantaneously captures full absorbance spectra from 220 - 1000 nm at resolutions of 1 - 10 nm.

Full absorbance spectra can be measured in less than one second per well - significantly faster than any conventional method. Alternatively, partial spectra or up to eight discrete wavelengths can be measured simultaneously.

The spectrometer incorporates a highly efficient optical grating and a solid state array detector. This allows the measurement of light intensity throughout the UV and visible parts of the spectrum. Similar to a monochromator, but much faster, the spectrometer can capture the entire UV/Vis spectrum of a sample within a second - no scanning needed.



From a single spectrum per well to spectra overlay plots

Standard Cuvette Port

More than just a microplate reader, the SPECTROstar^{Nano} has a standard cuvette port for single sample absorbance measurements. A range of sample sizes (0.3 μ L to 5 mL) can be easily measured with several special cuvettes. The SPECTROstar^{Nano} redefines cuvette-based kinetic assays. Full spectrum data can be collected in one second, allowing for polychromatic evaluation over the whole UV to visible spectrum.

The cuvette port door can be kept open during measurements, which is ideal for a flow-through cell. This advantage also allows readings to be taken before, during, and after manual addition of reagents.




Cuvette port for standard and micro cuvettes

LVis Plate

When quantitating DNA or RNA, using the least amount of sample is always preferable. With the LVis Plate, BMG LABTECH's low-volume microplate, sixteen individual samples of 2 μ L can be measured. The LVis Plate contains a tip rest for easy sample loading when using an eight channel pipettor. The LVis Plate also contains NIST traceable optical density filters and a holmium oxide filter for instrument quality control and performance testing.

In addition, a horizontal cuvette position can hold standard rectangular cuvettes for kinetic or endpoint studies.



LVis Plate for low-volume (2 μ L) measurements and instrument performance testing

Incubation, Shaking, and Gas Vent

The SPECTROstar ^{Nano} has incubation, multi-mode shaking capabilities, and a gas vent. Since all cell-based assays require the correct environment to be properly measured, the gas vent can be used to purge the microplate chamber with different gases. For example a mixture of O₂/CO₂ can be used to measure mammalian tissue cells. The built-in incubator uniformly heats the microplate chamber and the cuvette port up to 45°C, thereby allowing the optimal temperature to be used for all cellular growth assays.

The built-in incubator can also be used for protein denaturation studies by incrementally increasing the temperature followed by full spectrum readings. Multi-mode shaking - linear, orbital, and double orbital - permits flexibility when designing assays. Linear shaking allows for the fastest and most intense shaking conditions for cells that tend to clump. Double orbital shaking allows for the slowest and mildest conditions for sensitive cells. Orbital shaking is a cross in between the two.

Endpoint, Kinetics, and Well Scanning

Users can capture an endpoint signal that happens in seconds, or measure the kinetics of bacterial growth over time. In well scanning mode, the SPECTROstar ^{Nano} can easily deal with nonhomogeneous samples such as adherent cells by taking multiple measurements in each well with up to 900 data points/well. The MARS software displays each scan point graphically and creates a map for each well.

Automation


All BMG LABTECH microplate readers have similar dimensions and plate in/out locations. This minimizes the cost of automation solutions. Small footprint, automation-friendly plate carrier, and multiple robotic software interfaces allow for easy integration into all robotic platforms. Furthermore, having 1536-well capability allows for high-throughput full spectrum analysis.

Control and MARS Data Analysis Software

BMG LABTECH's software package provides an extensive range of possibilities for both test protocol definitions and data analysis. It is fully compliant with FDA regulation 21 CFR Part 11. The Control software defines instrument parameters and test protocols; the MARS software offers various tools for data analysis. Well-organized, versatile, easy to use, and powerful are just a few of the ways the MARS Data Analysis Software package is described by users. MARS provides several options to display data in a clear and concise format.

Data can be processed with powerful predefined templates or by using an extensive range of data calculation features. Examples include the automatic calculation of enzyme kinetic parameters (V_{max} and K_m) with a variety of fits based on Michaelis-Menten or Lineweaver-Burk equations, or the generation of standard curves based on the following curve fitting algorithms to calculate e.g. EC₅₀, IC₅₀, and r^2 values:

- Linear regression fit
- 4 parameter fit
- Point to point fit
- Segmental regression fit
- Cubic spline fit
- 2nd and 3rd polynomial fit
- Enzyme kinetic (e.g. Michaelis-Menten)



The MARS Data Analysis Software provides sophisticated tools for automated data reduction

Applications Center


Numerous publications for absorbance measurements such as application notes, posters, and scientific papers exemplify the versatility of the SPECTROstar^{Nano}.

A wide range of possible applications include:


- 260/280 ratio for DNA and RNA
- Protein quantitation
- ELISA
- Cell-based assays
- Enzyme activity assays

The versatility and flexibility of the SPECTROstar^{Nano} spectrometer are illustrated by the following examples:


- Kinetic studies on the metallation of porphyrin
- Protein determination with the Bradford assay
- Monitoring of bacterial growth at elevated hydrostatic pressure



Changes in visible spectrum accompanying zinc metallation of TPP. Arrows indicate the evolution of the absorption bands with time¹



3D well map of a Bradford absorbance spectrum taken over 60 seconds. The 475 nm peak decreases and the 595 nm peak increases



Dose dependant growth curves of *Candida albicans* in Tree Tea Oil (TTO). The mean of the quadruplicate data was taken and standard error bars presented.

BMG LABTECH continuously works with all major reagent companies to develop protocols and to optimize instrument settings for their existing assays and their newest kits. Visit the online Applications Center at www.bmglabtech.com to download all the leading applications, listed as:

- Application notes
- Scientific posters
- Peer-reviewed papers

The applications database reflects the expertise acquired by BMG LABTECH over the last 20 years. Over 3,800 published entries of peer-reviewed articles, application notes, and scientific posters demonstrate the flexibility and versatility of our readers, and their use in chemical and biological sciences.

Support and Training

BMG LABTECH operates globally through an extensive network of subsidiaries and well-trained distributors. Customers can rely on PhD-level support and assistance with regard to software, assay development, or general enquiries related to the SPECTROstar^{Nano} and all other BMG LABTECH microplate reading solutions.

^{1, 2} The graphs were taken from BMG LABTECH's Application Notes AN 178 and AN 191.

SPECTROstar^{Nano} - Technical Specifications



Due to the modularity of BMG LABTECH's instruments, all or combinations of the features below can be installed at purchase or upgraded at any time. Please contact your local representative for more details or a quote.

Detection Modes	UV/Vis absorbance spectrum
Measurement Modes	Endpoint and kinetic measurements Well scanning
Microplate Formats	6- to 1536-well plates, user-definable
Microplate Carrier	Robot compatible
Cuvette Port	Cuvette port for cuvettes with 10 mm path length Micro cuvettes, Traycell compatible Beam height 8.5 mm
Light Source	High-energy xenon flashlamp
Detector	Spectrometer with CCD
Spectral Range	220 - 1000 nm
Sensitivity	Selectable spectral resolution: 1, 2, 5, 10 nm OD range: 0 to 4 OD Accuracy: < 1% at 2 OD Precision: < 0.5% at 1 OD and < 0.8% at 2 OD
Read Times	Full spectrum from 220 to 1000 nm in less than 1 sec/well
Shaking	Linear, orbital, and double orbital User-definable time and speed
Incubation	+3°C above ambient up to 45°C for microplate and cuvette
Gas Vent	System to inject an atmosphere or to pull a vacuum into the reader
Software	Multi-user software package including Reader Control and MARS Data Analysis Software, FDA 21 CFR Part 11 compliant
Computer Interface	USB 2.0 compatible to USB 1.1
Dimensions	Width: 36 cm, depth: 50 cm, height: 16 cm Weight: 10 kg
Accessories	
LVis Plate	Microplate designed to measure 16 low-volume (2 µL) samples and standard cuvettes. Incorporates NIST-traceable filters and holmium oxide standards for instrument performance testing. Sensitivity: 2 ng/µL dsDNA
Barcode Reader	External manual barcode reader
THERMOstar	Microplate incubator and shaker

Specifications are subject to change.

Limit of detection (sensitivity) was calculated according to the IUPAC standard: $3 \times (SD_{blank}) / \text{slope}$
SPECTROstar is a registered trademark of BMG LABTECH GmbH.

© 2014 All rights reserved. All logos and trademarks are the property of BMG LABTECH.



Made in Germany

Headquarters Germany
BMG LABTECH GmbH
Allmendgrün 8
77799 Ortenberg
Tel. +49 781 96968 -0
germany@bmglabtech.com

Australia
BMG LABTECH Pty. Ltd.
2/24 Carbine Way
Mornington, Victoria, 3931
Tel. +61 3 5973 4744
australia@bmglabtech.com

France
BMG LABTECH SARL
7, Rue Roland Martin
94500 Champigny s/Marne
Tel. +33 1 48 86 20 20
france@bmglabtech.com

Japan
BMG LABTECH JAPAN Ltd.
1-6-2, Shimo-cho
Omiya-ku
330-0844 Saitama City
Tel. +81 48 647 7217
japan@bmglabtech.com

UK
BMG LABTECH Ltd.
5 Alton House Office Park
Gatehouse Way
Aylesbury
HP19 8YB
Tel. +44 1296 336650
uksales@bmglabtech.com

USA
BMG LABTECH Inc.
13000 Weston Parkway
Suite 109
Cary, NC 27513
Tel. +1 877 264 5227
usa@bmglabtech.com



Go Green Policy

The SPECTROstar^{Nano} power supply uses less than 0.5 watts of power when the instrument is in standby mode (Efficiency Level V). With EU RoHS construction standards this instrument complies with the EU Restriction of Hazardous Substances Directive (RoHS) and has the smallest impact on the environment.

www.bmglabtech.com