



RPFR05043 Solid Door Pharmacy & Vaccine Refrigerator*

| Part Number | RPFR05043 |
|-----------------------------------------------------------------|------------------------------------------------|
| Capacity (Litres) | 150 |
| Style | Under Counter |
| Exterior Dimensions (HxWxD) | 820 x 600 x 620mm |
| Interior Dimensions (HxWxD) | 680 x 540 x 350mm |
| Energy Consumption | 3.75 kWh/24h |
| Operating Temperature | +4°C |
| Temperature Set Point Range | +2°C to +8°C |
| Interior Construction | Stainless Steel |
| Exterior Construction | Stainless Steel |
| Shelves | 2 |
| Shelf dimensions (W x D) | 530 x 325mm |
| Door Hinging | Right hand only |
| Door Style | Solid (Glass Door Model available RPF05043) |
| Internal Light | Glass Door Model only |
| Rear Castors | ✓ |
| Door Lock | ✓ |
| Access Port | Optional Extra |
| Digital Temperature Control | ✓ |
| Digital Temperature Display | ✓ |
| High Temperature alarm | ✓ |
| Low Temperature alarm | ✓ |
| Door Open alarm | ✓ |
| Auto defrost | ✓ |
| Fan Assisted Cooling | ✓ |
| Environmentally friendly insulation | ✓ |
| Eco Friendly Refrigerant | ✓ |
| Ammonia Free | ✓ |
| DH Green Book, RPS/ GPhC, NaTHNaC and RCVS Guidelines compliant | ✓ |

Robustly constructed from Stainless Steel this refrigerator from the **Labcold** Advanced range is designed to provide close temperature control for temperature sensitive products.



The fan circulated interior has no internal 'step', making it easy to clean and stock and to allow for maximum air flow for even temperature distribution. The two shelves are fully adjustable, so it is possible to maximise the space available.

Insulated with Envirofoam, which is derived from plant rather than petrochemical oils, the refrigerant is also ammonia, HCFC and CFC free.

Fitted with a sturdy door lock, this Advanced refrigerator is ideal for busy environments such as busy laboratories or pharmacies.

* Glass door model available. Please order RPF05043.