

The SciPette Aspi-Vac Safety Vacuum Aspiration System is ideal for the safe removal of hazardous waste such as the removal of media from cell cultures, supernatants from centrifugation or wash solutions from immune-assay well-plates.

- ◆ Integrated pump with vacuum control for easy operation
- ◆ Self-closing connectors for quick and easy installation, plus the prevention of drips and aerosol escape
- ◆ Autoclavable waste collection bottle and nozzle
- ◆ Level sensor prevents overflow
- ◆ 4 Litre waste collection bottle

Order Number	Description
SP3080	SciPette Aspi-Vac Safety Vacuum Aspiration System

Bottle Volume	4L
Vacuum Range	0 to -600 mbar (adjustable)
Pump Flow Rate (rated value)	15 L/min (air)
Aspiration Rate of Liquids	17 ml/s (aspiration pipette)
Dimensions (WxDxH)	180 x 240 x 450 mm
Weight	3 kg
Input Voltage	100-240 V, 50/60Hz

