

The LETO is the SciPette TITAN's little sister, featuring the same great battery life, robust construction and universal nose cone, the LETO has been simplified and is more elegant. With enhanced ergonomics, the LETO is extremely comfortable for both left and right handed users.

In place of the TITAN's digital screen and 8 set speeds, the LETO features a continuously adjustable thumb wheel to set the ideal aspiration and dispensing speed.



- ◆ **Light weight electronic pipette filler**
- ◆ **One handed operation**
- ◆ **Accepts stranded pipettes from 0.1ml – 100ml**
- ◆ **Variable speed**
- ◆ **Lithium-ion battery offers long runtime per charge and no battery memory problems**
- ◆ **Hydrophobic filter 0.45µm**
- ◆ **Fully autoclavable nose cone**
- ◆ **Complete with charger, spare filter and holder**
- ◆ **Elegant design**

| Part Number | SP2222 |
|----------------------|---|
| Aspiration Speed | Variable, analogue thumb wheel |
| Dispense Speed | Variable, analogue thumb wheel |
| Battery | Lithium-ion |
| Battery Service Life | More than 8 hours of intermittent use |
| Charging Time | 2-3 hours |
| Pipette Types | Glass or Plastic pipette 0.1-100ml, Pasteur |
| Filter | 0.45 µm hydrophobic |