


SciPette

Liquid Handling Solutions

TITAN Pipette Filler



www.sciquip.co.uk



Fully autoclavable nosecone housing simplifies disinfecting and decontamination

Large LCD display provides visual confirmation of remaining battery charge and speed settings

Select the desired speed using the speed adjustment control

Separate aspirate and dispense speed controls

SciPette TITAN

- Full volume range 0.1 - 100ml
- Easy to use one-handed operation
- Lightweight, ergonomic design allows longer, fatigue-free pipetting
- Self-regulate the aspirate and dispense speed by circuit feedback
- Lithium-ion battery offers long runtime per charge and no battery memory problems
- Low battery indication
- Fills 100ml quickly with powerful motor

Technical Specifications

Part No.	SP2200
Aspirate Speeds	8
Dispense Speeds	8 + Gravity Dispense
Battery	Lithium-ion
Battery Service Life	More than 15 hours of continuous use
Charging Time	Approx 2hrs
Pipette Types	Glass or Plastic Pipettes (0.1-100ml), Pasteur pipettes
Filter	Hydrophobic

SciQuip Ltd
Newtown
Wem
Shrewsbury
Shropshire
SY4 3QJ, UK

Tel: +44 (0)1939 234222
Fax: +44 (0)1939 234221
Email: info@sciquip.co.uk
www.sciquip.co.uk