

SciPette Pipette Pumps

SciPette Pipette Pumps are a low cost, syringe-style pipetting device, that offers one-hand control of filling and delivery for pipettes that are calibrated to the tip. The feature of the side lever enables complete liquid dispensing, without a droplet remaining.

Designed for single-handed operation, The Pipette Pump identifies itself with a thumbwheel, allowing precise aspiration or dispensing. The controller can be easily disassembled for cleaning.

Technical Specifications:

- Comfortable design
- Precise pipetting control
- Rapid delivery lever
- Easily disassembled for cleaning
- Color coded for maximum volumes
- Volume capacity: 2ml, 10ml and 25ml

