

SciQuip 3D Shakers

The SciQuip laboratory benchtop 3D shakers create a perfectly smooth circular tilting motion within culture flasks and petri dishes. Both units are fitted with quiet, powerful and maintenance free DC motors and are a reliable and robust solution for any laboratory. The range covers the following key features:

- Run Timer (Pro only)
- Soft start with smooth acceleration
- Wide speed range for a variety of applications
- Supplied with grip mat platform
- Cold room and incubator compatible
- IP21 rated

Technical Specifications

	3D (Basic)	3D (Pro)
Order Number	SP2258	SP2254
Shaking movement	3D Orbiting	3D Orbiting
Tilting angle		
Max. load	3.0kg	5.0kg
Speed range	0-80	10-70
Display	No	LCD
Control	Analogue	Digital
Timer	No	Yes
Time setting range	---	1min-20h
Dimensions (WxHxD)	270x130x330mm	360x106x430mm
Weight	2.7	9
Power	20W	50W
Platform Fitted as standard	Universal Grip Mat	Universal Grip Mat
Other platforms available	Yes	Yes