

GyroStir LCD Digital Magnetic Stirrer SP2255-MS6

- ◆ Max. heating temperature is 380°C
- ◆ High resolution LCD displays actual temperature and speed
- ◆ Brushless DC motor is maintenance free
- ◆ Aluminum cover with ceramic work plate - immediate heat transfer
- ◆ External temperature control is possible with temperature sensor
- ◆ PT1000
- ◆ Digital speed control with max. speed up to 1500rpm
- ◆ Max. stirring quantity of H₂O at 5L
- ◆ 2 Year Warranty



GyroStir Glass Ceramic Square Magnetic Stirrer SP2245

- ◆ 7" ceramic acid resistant work plate
- ◆ 0-1500rpm speed range
- ◆ Stirrers up to 10L of water
- ◆ Powerful magnet
- ◆ Simple analogue control
- ◆ 2 year warranty



GyroStir Basic Digital Magnetic Stirrer

- ◆ Digital Magnetic Stirrer
- ◆ ABS acid resistant* plastic 5" plate
- ◆ Speed range 100 – 1,500 rpm
- ◆ Max stirring capacity 3 liters (water)
- ◆ LED display – both set & real time values
- ◆ Single one-knob control
- ◆ 2 year warranty

*weak acids and alkalies only



GyroStir 10 Multi Position Magnetic Stirrer

- ◆ Anti-slip chemical resistant surface
- ◆ High-performance multi-position magnetic stirrer with 10 stirring position, optimum use of laboratory space
- ◆ Stainless steel surface allowing easy cleaning and providing protection against penetration of liquids
- ◆ Simultaneously operating stirrers, absolute consistency over sample condition for the individual samples
- ◆ Stepless speed regulation. Wide speed range: 0-1100 rpm/min
- ◆ Motor with electronic speed control, constant speed even during changes in load
- ◆ 2 year warranty



SciQuip Ltd
Newtown
Wem
Shropshire
SY4 5NU, UK

Tel: +44 (0)1939 234222
Fax: +44 (0)1939 234221
Email: info@sciquip.co.uk
www.sciquip.co.uk

Ordering Information

Part Number	Description
-------------	-------------

SP2236-FS	GyroStir FlatSpin 3.5" Magnetic Stirrer
-----------	---

SP2201-B-A	GyroStir Basic – Analogue 5" Magnetic Stirrer
------------	---

SP2201-B-D	GyroStir Basic – Digital 5" Magnetic Stirrer
------------	--

SP2240-10A	GyroStir 10 Position Magnetic Stirrer
------------	---------------------------------------

Stirrer Accessories

Part Number	Description
-------------	-------------

SP-MSBAR10	Magnetic Stirrer bar (10mm x 6mm)
------------	-----------------------------------

SP-MSBAR15	Magnetic Stirrer bar (15mm x 8mm)
------------	-----------------------------------

SP-MSBAR20	Magnetic Stirrer bar (20mm x 8mm)
------------	-----------------------------------

SP-MSBAR25	Magnetic Stirrer bar (25mm x 8mm)
------------	-----------------------------------

SP-MSBAR30	Magnetic Stirrer bar (30mm x 6mm)
------------	-----------------------------------

SP-MSBAR40	Magnetic Stirrer bar (40mm x 8mm)
------------	-----------------------------------

SP-MSBAR50	Magnetic Stirrer bar (50mm x 8mm)
------------	-----------------------------------

SP-MSBAR65	Magnetic Stirrer bar (65mm x 8mm)
------------	-----------------------------------

SP-MSBAR80	Magnetic Stirrer bar (80mm x 13mm)
------------	------------------------------------

SP-MARRM	Magnetic Stirrer bar remover, length of 200mm
----------	---

Model	GyroStir FlatSpin	GyroStir Basic Analogue	GyroStir Basic Digital	GyroStir 10
Part number	SP2236-FS	SP2201-B-A	SP2201-B-D	SP2240-10A
Number of positions	1	1	1	10
Plate size	90mm	135mm	135mm	180x450mm
Max H2O	0.8L	3L	3L	0.4L per position
Plate material	Plastic	Plastic	Plastic	Stainless Steel & Silicone film
Display	None	None	LCD	None
Max Length Stirrer Bar	30mm	50mm	50mm	40mm
Dimensions DxWxH	120x175x15mm	150x260x80mm	150x260x80mm	182x622x65mm
Ingress protection rating	IP65	IP42	IP42	IP42

The GyroStir Cell Culture Magnetic Stirrer is a heady duty work horse specifically designed for laboratories processing cell cultures. Available in both single and 4 position, the platforms are large enough for 3L spinner and water jacketed culture flasks and their intuitive programmability and all-steel construction make them the ideal choice for any lab.

- Specifically designed for low speed routine process cell cultures
- Intuitive controller
- Large, clear LED display
- Many program modes such as continuous run, cycling, alternate speed etc.
- Robust steel construction
- Single and 4 position versions
- 2 year warranty



	GyroStir Cell Culture – 1 Position	GyroStir Cell Culture – 4 Position
Order Number	SP2238-1	SP2238-4
Speed Range	5 – 200 rpm	5 – 200 rpm
Max Bottle Volume	3L (x1)	3L (x4)
Power	15W	15W
Voltage	100-240V (50-60Hz)	100-240V (50-60Hz)
Dimensions (WDH)	205 x 250 x 90mm	390 x 450 x 90mm
Warranty	2 Years	2 Years

