

SciQuip Magnetic Stirrers up to 3L

DL2201-B-DA

LED Digital Magnetic Stirrer

Features

- Digital speed control within the range of 0-1500rpm
- Max. stirring quantity of H₂O at 3L
- LED display shows speed
- Nylon+GF housing provides a good performance of chemical resistance



DL2201-B-A

Magnetic Stirrer

Features

- A wide range of 0-1500rpm
- Max. stirring quantity of H₂O at 3L
- Nylon+GF housing provides a good performance of chemical resistance

Specifications

	DL2201-B-DA	DL2201-B-A
Work plate dimension	φ135mm(5 inch)	φ135mm(5 inch)
Work plate material	Nylon+GF	Nylon+GF
Motor type	Brushless DC motor	Brushless DC motor
Stirring positions	1	1
Max. stirring quantity[H ₂ O]	3L	3L
Max. magnetic bar[length]	50mm	50mm
Speed range	0-1500rpm	0-1500rpm
Speed display	LED	Scale
Protection class	IP42	IP42
Motor input power	5W	5W
Motor output power	3W	3W
Power	15W	10W
Voltage	100-120/200-240V,50/60Hz	100-120/200-240V,50/60Hz
Dimension[WxDxH]	150x260x80mm	150x260x80 mm
Weight	1.8kg	1.8kg
Permissible ambient temperature and humidity	10-40°C, 80%RH	10-40°C, 80%RH

DL2236-FS

Ultra Flat Compact Magnetic Stirrer

With excellent stirring functionality, the company shape of the FlatSpin takes less bench space and makes storage easier

A wide range of speed from 15 to 1500rpm
Automatic reverse rotation for excellent mixing
Protection class (IP65)
PET top provides excellent chemical resistance
Ultra silence and vibration-free



DL2236-ECO

Economical Magnetic Stirrer

EcoStir is a magnetic stirrer customized for liquid mixing, suitable for chemical analysis, biology reagent mixing, regular laboratory experiments and liquid handling.

Maintenance free brushless DC motor
Rare earth magnets give powerful stirring
Stepless speed regulation 300-2000rpm
Ultra silence and vibration-free design.
PET top provides excellent chemical resistance
Shaft can be fitted with a PH sensor or sensors for other purposes



Specifications	DL2236-FS	DL2236-ECO
Work plate Dimension	φ90mm	φ120mm
Work plate material	PET	PET
Motor type	Motorless	Brushless DC motor
Stirring positions	1	1
Max. stirring quantity[H ₂ O]	0.8L	1.5L
Max. magnetic bar[length]	40mm	40mm
Speed range	15-1500rpm	300-2000rpm
Speed setting	button	scale
Protection class	IP65	IP21
Motor output power	3W	1.2W
Power	5W	7.2W
Voltage	100-120/200-240V,50/60Hz	100-120/200-240V,50/60Hz
Dimension [W x D x H]	175x120x15mm	150x145x50mm
Weight	0.3kg	0.3kg
Permissible ambient temperature and humidity	5-40°C,80%RH	5-40°C,80%RH

SciQuip Ltd
Newtown
Wem
Shrewsbury
Shropshire
SY4 5NU, UK

Tel: +44 (0)1939 234222
Fax: +44 (0)1939 234221
Email: info@sciquip.co.uk
www.sciquip.co.uk