

- ◆ 6 roller bars – 5 position
- ◆ Bars roll and rock with 24mm off-center cams at the rear
- ◆ Gentle and highly effective agitation for preventing blood coagulation and maintaining liquid-solid suspensions
- ◆ Wide speed range for all application
- ◆ High torque maintenance free DC motor
- ◆ Suitable for use in incubators and cold rooms
- ◆ 2 year warranty



The SciQuip Tube Roller Pro also includes

- ◆ Timed run function
- ◆ LCD display
- ◆ Full Digital control



	SciQuip Tube Roller Basic	SciQuip Tube Roller Pro
Part Number	SP2270	SP2272
Number of rollers	6	6
Roller length	280mm	280mm
Load Capacity	4kg	4kg
Motion	Rocking & Rolling	Rocking & Rolling
Rock Amplitude	24mm	24mm
Speed Range	0-70 rpm	10-70 rpm
Timer function	No	Yes
Display	No	Yes
Operation Mode	Continuous	Continuous / Timed
Ingress Protection	IP 21 DIN EN60529	IP 21 DIN EN60529
Dimensions (DWH)	450 x 260 x 120mm	450 x 260 x 120mm
Power	25W	30W