


- SciQuip Disc Rotators are supplied with a disc as standard
- Smooth motion provides gentle and effective tube mixing
- Large speed range
- Ideal for microtubes, 15ml and 50ml conical culture tubes
- Compact design
- High torque, zero maintenance DC motor
- Suitable for use in cold rooms and incubators
- 2 year warranty

The SciQuip Disc Rotator Pro also features:

- Twin disc capability offering double the capacity of its competitors
- LCD display
- Timer option
- Adjustable angle
- Fully digital control



Technical Specification

	Disc Rotator Basic	Disc Rotator Pro
Part Number	SP2278	SP2274
Standard Equipment	Disc & 12 of each clamp size	Disc with 8 x 50ml clamps
Max. Number of Discs	1	2
Max. Capacity	48 x 1.5/2ml 16 x 15ml 8 x 50ml	60 x 1.5/2ml Per Disc 16 x 15ml Per Disc 8 x 50ml Per Disc
Speed Range	0-80 RPM	10-70 RPM
Disc Angle	90° Fixed	0-90° Variable
Controller Mode	Continuous Run	Continuous or Timed Run
Ingress Protection	IP 21 DIN EN60529	IP 21 DIN EN60529
Dimensions (DWH)	300x220x310mm	280x210x300mm
Power	20W	40W