


SciQuip Overhead Tube Rotator

- Supplied with rotisserie adapter as standard
- Ideal for 1.5/2.0ml microtubes, 15ml and 50ml conical culture tubes
- Large speed range
- Compact design
- High torque, zero maintenance DC motor
- Suitable for use in cold rooms and incubators
- 2 year warranty



The SciQuip Overhead Rotator Pro also features:

- LCD display
- Timer option
- Adjustable angle
- Fully digital control



	Disc Overhead Rotator Basic	Disc Overhead Rotator Pro
Part Number	SP2279	SP2276
Standard Equipment	Rotor & 12 of each clamp size	Rotor with 16 x 50ml clamps
Max. Capacity	48 x 1.5/2ml 24 x 15ml 12 x 50ml	48 x 1.5/2ml 24 x 15ml 24 x 50ml
Speed Range	0-80 RPM	10-70 RPM
Disc Angle	Fixed	0-90° Variable
Controller Mode	Continuous Run	Continuous or Timed Run
Ingress Protection	IP 21 DIN EN60529	IP 21 DIN EN60529
Dimensions (DWH)	150x530x190mm	220x510x260mm
Power	20W	40W

SciQuip Ltd
Newtown
Wem
Shropshire
SY4 5NU, UK

Tel: +44 (0)1939 234222
Fax: +44 (0)1939 234221
Email: info@sciquip.co.uk
www.sciquip.co.uk