

SciQuip Incu C-Series

The SciQuip Incu C-series cooled incubators incorporate a precise cooling system providing excellent temperature uniformity throughout making them ideal for any application requiring accurate below ambient temperature control. All parameters are displayed in the large LCD screen for fast easy operation and read-out together with a graphic program display. The SciQuip Incubators can store up to 9 programs of 18 steps, each program can be made to repeat from 1 to 99 times automatically.

- Capacity 90ltr – 430ltr
- Temp 4°C to 65°C
- Fully programmable
- Electro-polished chamber & shelves
- Refrigerated models come with a side access port
- Refrigerated models come with a internal lighting
- Audible & visual alarms
- Password protection
- Inner glass door
- Window in door as standard
- 50mm port as standard
- Interior lamp as standard
- **All units are supplied with 2 shelves, as standard**



Model	SQ-4650	SQ-4655	SQ-4660	SQ-4665
Volume (L)	90	160	270	430
Temperature Range	Ambient +4°C to 65°C	Ambient +4°C to 65°C	Ambient +4°C to 65°C	Ambient +4°C to 65°C
Temperature Accuracy	±0.01°C	±0.01°C	±0.01°C	±0.01°C
Temperature Uniformity	±1.0°C @ 37°C	±1.0°C @ 37°C	±1.0°C @ 37°C	±1.0°C @ 37°C
Alarm	Enabled	Enabled	Enabled	Enabled
Timer	0-999 minutes	0-999 minutes	0-999 minutes	0-999 minutes
Setting	Digital	Digital	Digital	Digital
Display	LCD	LCD	LCD	LCD
Shelf Included	2 (max 11)	2 (max 15)	2 (max 18)	2 (max 25)
Shelf Size (mm) (WxD)	310x356	410x456	513x556	555x656
Inner Dimensions (mm) (WxD)	400x410x550	500x500x650	600x600x750	700x645x950
Exterior Dimensions (mm) (WxDxH)	550x555x1280	630x740x1380	750x840x1640	910x950x1840
Packing Dimensions (mm) (WxDxH)	620x625x1440	700x810x1540	820x810x1640	910x950x1840
Net/Gross Weight (KG)	68/108	98/145	130/180	180/220
Power (W)	710	860	950	1350
Electricity	220-240V 50/60 Hz	220-240V 50/60 Hz	220-240V 50/60 Hz	220-240V 50/60 Hz
Approval	CE, ISO	CE, ISO	CE, ISO	CE, ISO
Warranty	2 Years	2 Years	2 Years	2 Years

SciQuip Ltd
Newtown
Wem
Shrewsbury
Shropshire
SY4 5NU, UK

Tel: +44 (0)1939 234222
Fax: +44 (0)1939 234221
Email: info@sciquip.co.uk
www.sciquip.co.uk