

SciQuip Overstock Items - KERN

WTB 1K-4N Kern Bench Scale 00002kg ; 15kg



With an innovated coloured changing display and IP65 rated this compact set of scales is ideal little balance for the food industry!

- Large backlit LCD display, digit height 25 mm
- Weighing plate dimensions, stainless steel, W×D×H 260×200×16,2 mm
- Overall dimensions W×D×H 285×315×128 mm
- Net weight approx. 3,2 kg

SD-EOB 35K10 Kern Bench Scale 10g ; 35kg



These scales are ideal to be used for parcels and alternative a great tool for veterinarians. Extremely easy to use, with 4 key operation and ensure accuracy with every measurement by automatically filtering out vibrations so perfect for animals that may be a little restless!

- Large LCD display, digit height 25 mm
- Dimensions of weighing plate W×D×H 315×305×65 mm
- Cable length of display device approx. 1,8 m
- Optional battery operation, 6 x 1.5 V AA not included, operating time up to 60 h

SD-MPD 250K100M Kern Personal floor scale 0,1kg ; 250kg



Ideal for doctor's surgeries, chemists and hospitals: easy to use and robust with a non-slip and wear resistant surface.

- Approved for medical use according to 93/42/EEC
- The bestseller for doctor's surgeries, chemists and hospitals: easy to use with no time wasted
- Simply tap on it, stand on it, read the weight
- For easy and hygienic cleaning

SD-MPC 250K100M Kern Personal floor scale Max 250kg; e=0,1kg ; d=0,1kg



Ideal for doctor's surgeries, chemists and hospitals: easy to use and robust with a non-slip and wear resistant surface. Additional Mother and child function.

- Approved for medical use according to 93/42/EEC
- The bestseller for doctor's surgeries, chemists and hospitals: easy to use with no time wasted
- Simply tap on it, stand on it, read the weight
- For easy and hygienic cleaning

SciQuip Overstock Item

SD-IFB 6K-4A Kern Industrial platform scale Max 6000g; d=200mg



Kerns Platform scales suitable for tougher industry use, with a stainless-steel plate and IP65 rated.

- Large backlit LCD display, digit height 52 mm
- Weighing plate dimensions, stainless steel W×D×H 300×240×110 mm
- Cable length of display device approx. 3 m
- IP rated against dust and water splashes

SD-ABT 220-4NM Kern Analytical Balance ABT. Max 220g x 0.1 readability



The Kern premium model with single-cell weighing system and internal calibration, this is a great addition and good price!

- Large LCD display, digit height 14 mm
- Dimensions weighing surface, stainless steel, Ø 80 mm
- Overall dimensions (incl. draught shield) W×D×H 217×356×338 mm
- Weighing space W×D×H 168×172×223 mm

SD-EW 820-2NM Kern Precision Lab Balance 820g x 0.01g



This classic balance with robust tuning fork measuring system

- Large LCD display, digit height 17 mm
- Weighing plate dimensions, stainless steel Ø 118 mm WxD 170x142 mm
- Overall dimensions (without draught shield) WxDxH 182x235x75 mm
- Externally calibrated

SD-ADB600-C3 Kern Analytical balance Max 120 g; d=0,1 mg



The high level of accuracy saves hard cash wherever you are weighing valuable precious stones

- Adjusting program CAL for quick setting of the balance accuracy using an external test weight
- Level indicator and levelling feet for precise levelling of the scale, fitted as standard, to give the most accurate weighing result
- Large glass draught shield with 3 sliding doors for easy access to the items being weighed
- Compact size, practical for small spaces

SciQuip Overstock Items

SD-EMS12K1 Kern EMS Precision balance 1g ; 12000g



Entry level model in the low-cost range with large weighing plate

- Especially suitable for use in schools and universities, for example for biology, chemistry, physics
- Large, shock proof weighing plate made of plastic
- Particularly flat design
- Ergonomically optimised key pad with large keys and a high-contrast LCD display
- Secure and non-slip positioning with rubber feet

SD-HCB50K100 Kern HCB Hanging Scales 100 g ; 50 kg



The handy scales for higher loads up to 1 t

- Ideal for rapid control in goods-in and goods-out
- Essential too in the private sector to determine the weight of fish, game, fruits, bicycle parts, suitcases etc.
- Hold function: For easy reading of the weighing result, the display can be "frozen" in different ways. Either automatically when the weighing value remains unchanged or manually by pressing the Hold key
- Peak load display (peak hold), measuring frequency 5 Hz

SD-HCB200K500 Kern HCB Hanging scale 0,5 kg ; 200 kg



The handy scales for higher loads up to 1 t

- Ideal for rapid control in goods-in and goods-out
- Essential too in the private sector to determine the weight of fish, game, fruits, bicycle parts, suitcases etc.
- Hold function: For easy reading of the weighing result, the display can be "frozen" in different ways. Either automatically when the weighing value remains unchanged or manually by pressing the Hold key
- Peak load display (peak hold), measuring frequency 5 Hz

SD-ADJ600-C3 Kern Analytical balance Max 120 g; d=0,1 mg



The price leader in analytical balances - now also with internal adjustment

- Adjusting program CAL for quick setting of the balance accuracy using an external test weight
- Level indicator and levelling feet for precise levelling of the scale, fitted as standard, to give the most accurate weighing result
- Large glass draught shield with 3 sliding doors for easy access to the items being weighed
- Compact size, practical for small spaces

SciQuip Overstock Items

SD-GAB 12K0.1N KERN GAB BENCH SCALES



Checkweighing and portioning scale, also with EC type approval [M]

- Large backlit LCD display, digit height 24 mm
- Dimensions weighing surface, stainless steel, W×D 294×225 mm
- Overall dimensions W×D×H 320×330×125 mm
- Net weight approx. 3 kg
- Permissible ambient temperature GAB-N: 0°C/40°C GAB-M: -10°C/40°C

SD-DE 120K10A Kern DE Platform scale 10 g ; 120 kg



The long lasting story of success with dust and spray protected display device

- High mobility
- Flexible positioning of display device
- Display device IP65
- Stainless steel weighing plate
- PRE-TARE function



The Scientific Equipment Specialists

Sustainability in Science

SciQuip Overstock Items



The Scientific Equipment Specialists

Sustainability in Science

SciQuip Overstock Items