

## SciQuip Overstock Items



### WTB 1K-4N Kern Bench Scale 00002kg ; 15kg

With an innovated coloured changing display and IP65 rated this compact set of scales is ideal little balance for the food industry!

- Large backlit LCD display, digit height 25 mm
- Weighing plate dimensions, stainless steel, WxDxH 260x200x16,2 mm
- Overall dimensions WxDxH 285x315x128 mm
- Net weight approx. 3,2 kg



### SD-EOB 35K10 Kern Bench Scale 10g ; 35kg

These scales are ideal to be used for parcels and alternative a great tool for veterinarians. Extremely easy to use, with 4 key operation and ensure accuracy with every measurement by automatically filtering out vibrations so perfect for animals that may be a little restless!

- Large LCD display, digit height 25 mm
- Dimensions of weighing plate WxDxH 315x305x65 mm
- Cable length of display device approx. 1,8 m
- Optional battery operation, 6 x 1.5 V AA not included, operating time up to 60 h



### SD-MPD 250K100M Kern Personal floor scale 0,1kg ; 250kg

Ideal for doctor's surgeries, chemists and hospitals: easy to use and robust with a non-slip and wear resistant surface.

- Approved for medical use according to 93/42/EEC
- The bestseller for doctor's surgeries, chemists and hospitals: easy to use with no time wasted Simply tap on it, stand on it, read the weight
- For easy and hygienic cleaning



### SD-MPC 250K100M Kern Personal floor scale Max 250kg; e=0,1kg ; d=0,1kg

Ideal for doctor's surgeries, chemists and hospitals: easy to use and robust with a non-slip and wear resistant surface. Additional Mother and child function.

- Approved for medical use according to 93/42/EEC
- The bestseller for doctor's surgeries, chemists and hospitals: easy to use with no time wasted
- Simply tap on it, stand on it, read the weight
- For easy and hygienic cleaning


## SciQuip Overstock Items



### SD-IFB 6K-4A Kern Industrial platform scale Max 6000g; d=200mg

Kerns Platform scales suitable for tougher industry use, with a stainless-steel plate and IP65 rated.

- Large backlit LCD display, digit height 52 mm
- Weighing plate dimensions, stainless steel WxDxH 300x240x110 mm Cable length of display device approx. 3 m
- IP rated against dust and water splashes



### SD-ABT 220-4NM Kern Analytical Balance ABT. Max 220g x 0.1 readability

These scales are ideal to be used for parcels and alternative a great tool for veterinarians. Extremely easy to use, with 4 key operation and ensure accuracy with every measurement by automatically filtering out vibrations so perfect for animals that may be a little restless!


- Large LCD display, digit height 25 mm
- Dimensions of weighing plate WxDxH 315x305x65 mm
- Cable length of display device approx. 1,8 m
- Optional battery operation, 6 x 1.5 V AA not included, operating time up to 60 h



### SD-EW 820-2NM Kern Precision Lab Balance 820g x 0.01g

Ideal for doctor's surgeries, chemists and hospitals: easy to use and robust with a non-slip and wear resistant surface.

- Large LCD display, digit height 17 mm
- Weighing plate dimensions, stainless steel Ø 118 mm WxD 170x142 mm
- Overall dimensions (without draught shield) WxDxH 182x235x75 mm
- Externally calibrated



### SD-ADB600-C3 Kern Analytical balance Max 120 g; d=0,1 mg

Ideal for doctor's surgeries, chemists and hospitals: easy to use and robust with a non-slip and wear resistant surface. Additional Mother and child function.

- Adjusting program CAL for quick setting of the balance accuracy using an external test weight Level indicator and levelling feet for precise levelling of the scale, fitted as standard, to give the most accurate weighing result
- Large glass draught shield with 3 sliding doors for easy access to the items being weighed Compact size, practical for small spaces

## SciQuip Overstock Items



### SD-EMS12K1 Kern EMS Precision balance 1g ; 12000g

Entry level model in the low-cost range with large weighing plate

- Especially suitable for use in schools and universities, for example for biology, chemistry, physics
- Large, shock proof weighing plate made of plastic
- Particularly flat design
- Ergonomically optimised key pad with large keys and a high-contrast LCD display
- Secure and non-slip positioning with rubber feet



### SD-HCB50K100 Kern HCB Hanging Scales 100 g ; 50 kg

The handy scales for higher loads up to 1 t

- Ideal for rapid control in goods-in and goods-out
- Essential too in the private sector to determine the weight of fish, game, fruits, bicycle parts, suitcases etc.
- Hold function: For easy reading of the weighing result, the display can be "frozen" in different ways. Either automatically when the weighing value remains unchanged or manually by pressing the Hold key
- Peak load display (peak hold), measuring frequency 5 Hz



### SD-HCB200K500 Kern HCB Hanging scale 0,5 kg ; 200 kg

The handy scales for higher loads up to 1 t

- Ideal for rapid control in goods-in and goods-out
- Essential too in the private sector to determine the weight of fish, game, fruits, bicycle parts, suitcases etc.
- Hold function: For easy reading of the weighing result, the display can be "frozen" in different ways. Either automatically when the weighing value remains unchanged or manually by pressing the Hold key
- Peak load display (peak hold), measuring frequency 5 Hz



### SD-ADJ600-C3 Kern Analytical balance Max 120 g; d=0,1 mg

The price leader in analytical balances - now also with internal adjustment

- Adjusting program CAL for quick setting of the balance accuracy using an external test weight Level indicator and levelling feet for precise levelling of the scale, fitted as standard, to give the most accurate weighing result
- Large glass draught shield with 3 sliding doors for easy access to the items being weighed Compact size, practical for small spaces



## SciQuip Overstock Items



### SD-GAB 12K0.1N KERN GAB BENCH SCALES

Checkweighing and portioning scale, also with EC type approval [M]

- Large backlit LCD display, digit height 24 mm
- Dimensions weighing surface, stainless steel, WxD 294x225 mm
- Overall dimensions WxDxH 320x330x125 mm
- Net weight approx. 3 kg
- Permissible ambient temperature GAB-N: 0°C/40°C GAB-M: -10°C/40°C



### SD-DE 120K10A Kern DE Platform scale 10 g ; 120 kg

The long lasting story of success with dust and spray protected display device

- High mobility
- Flexible positioning of display device
- Display device IP65
- Stainless steel weighing plate PRE-TARE function



### SD-DE 120K10A Kern DE Platform scale 10 g ; 120 kg

The long lasting story of success with dust and spray protected display device

- High Mobility
- Flexible positioning of display device
- Display device IP65
- Stainless steel weighing plate
- PRE-TAR function



### SD-EMB-200-2 School balance

SD-EMB-200-2 School balance

- Simple and convenient 2-key operation, making them ideal for use in schools and universities
- Tare function facilitates formulation work
- Particularly flat design
- Secure and non-slip positioning with rubber feet
- Ready for use: Batteries included. AUTO-OFF function to preserve the batteries (can be switched off)
- Ring-shaped draught shield standard, only for models with weighing plate size A, weighing space ØxH 96x35 mm
- Hook for underfloor weighing standard

## SciQuip Overstock Items



### SD-EMB-2000-0S Set School Balance

The space miracle – stackable school balance for experiment-based teaching

- Simple and convenient 2-key operation
- Tare function facilitates formulation work
- Particularly flat design
- Secure and non-slip positioning with rubber feet
- Stackable for space-saving storage
- A practical transport safety lock prevents the weighing technology from being overloaded when stacked, stored or during transportation



### SD-EMS 6K0.1 Precision balance

Entry level model in the low-cost range with large weighing plate

- Especially suitable for use in schools and universities, for example for biology, chemistry, physics
- Large, shock proof weighing plate made of plastic
- Particularly flat design
- Ergonomically optimised key pad with large keys and a high-contrast LCD display
- Secure and non-slip positioning with rubber feet
- Adjusting program CAL for quick setting of the balance accuracy, external test weights at an additional price, see test weights



### SD-EMS 3000-2 Kern EMS School balance 0,01 g ; 3000 g - NEW

Entry level model in the low-cost range with large weighing plate

- Especially suitable for use in schools and universities, for example for biology, chemistry, physics
- Large, shock proof weighing plate made of plastic
- Particularly flat design
- Ergonomically optimised key pad with large keys and a high-contrast LCD display
- Secure and non-slip positioning with rubber feet
- Adjusting program CAL for quick setting of the balance accuracy, external test weights at an additional price, see test weights



### SD-EWJ 6000-1M Precision balance (SG)

High-quality precision balance with automatic internal adjustment, verification optional

- Automatic internal adjustment, time-controlled every 2 h, guarantees high degree of accuracy and makes the balance independent of its location
- Capacity display: A bargraph display lights up to show how much of the weighing range is still available
- USB data interface for transferring weighing data to the PC, printer etc.
- Protective working cover included with delivery

## SciQuip Overstock Items



### SD-GAB 12K0.1N Kern Bench scale 0,1 g ; 12 kg

Checkweighing and portioning scale


- Compact size, practical for small spaces
- High mobility: thanks to rechargeable battery operation (optional), compact, lightweight construction, it is suitable for the use in several locations (production, warehouse, dispatch department, etc.)
- Weighing with tolerance range (checkweighing): a visual and audible signal helps with portioning, dispensing or grading
- Summation of weight values
- Protective working cover included with delivery



### SD-PCD-3000-2-1 Precision balance (SG)

High-resolution precision balance with removable display for maximum flexibility

- Laboratory balance with separate platform: Ideal when working in a glove bag or fume cupboard. Particularly practical for weighing toxic, volatile or contaminated substances
- PRE-TARE function for manual subtraction of a known container weight, useful for checking fill-levels
- Freely programmable weighing unit, e.g. display direct in special units such as length of thread g/m, paper weight g/m<sup>2</sup>, or similar
- Level indicator and levelling feet for precise levelling of the scale, fitted as standard, to give the most accurate weighing result
- Protective working cover included with delivery



### SD-PFB-200-3 Platform scale

Quick-display precision balance with user-friendly concept of operation mode – now with larger housing for more stability

- Compact size, practical for small spaces
- Capacity display: A bargraph display lights up to show how much of the weighing range is still available
- Level indicator and levelling feet for precise levelling of the scale, fitted as standard, to give the most accurate weighing result
- Draught shield standard, weighing space W×D×H 132×1125×65 mm