





Modification description for:

(SBS48)

XLX2000 G2

=>

XLX2000 G3

LVL'D UP
The next Generation



XLX2000 G2 (SBS48) => XLX2000 G3

Overall height identical

Working volume identical

Volume bottom point identical

XLX2000 G2 (SBS48) => XLX2000 G3




Tube height without cap at G3 is 0,5mm smaller than at G2

➔ No problem for switching from G2 to G3

XLX2000 G2 (SBS48) => XLX2000 G3

No dots on top of cap

Rounded ends of black component




Torsional securit 4x 90° with 4 in black area

2 dots on the top of the cap

Engraved logos on tube

Sharp corner ends of black component

Torsional securit 4x 90° with 2 in transparent area



XLX2000 G2 (SBS48) => XLX2000 G3

Both generations are compatible with the new automatic friendly generation of the rack



LVL'D UP
the next Generation

Contact

Waldemar Janzen
waldemar.janzen@lvl-technologies.com
Mob: +49 (0) 173 741 0392

OR

LVL technologies GmbH & Co. KG
info@lvl-technologies.com
Tel: +49 7951 95613-20

LVL'D UP
The next Generation