

High Speed Palm Micro Centrifuge

This centrifuges is ideal for quick spins, micro-filtration of samples, cell separations and many other routine laboratory procedures.

Top Features

- **Max Speed: 10,000 rpm (7490xg)**
- **OLED Display:** Easily switch speed and time settings
- **DC Motor Drive:** Ensures smooth and quiet operation
- **Safe and Durable Rotor:** Features clamp locking design for easy operation
- **Wide Voltage Technology:** Speed remains unaffected by voltage fluctuations, ensuring high-speed accuracy
- **Tube Compatibility:** Accommodates 0.5/1.5/2mL x 12 tubes and 4 x PCR8 strips, as well as 0.2mL x 44 PCR tubes
- **Screen Buttons:** Adjust motor start/stop with options for both soft braking and instant stop
- **Simple Design:** Lightweight at only 1.5 kg



Technical Specifications

Max Speed	10000rpm
Max RCF	6708xg-1.5/2mL; 6596-7490xg-PCR8 strips;
Rotor Capacity	0.2mLx44 PCR tubes 0.5/1.5/2mLx12+ 0.2mLPCR 8 stripsx4
Run Time	10s-20min
Display	OLED
Acceleration/Deceleration	1 ↑ / 2 ↓
Driving Motor	DC motor
Power	AC100-240V/50Hz/60Hz 40W
Noise Level	≤58dB
Dimension (D x W x H)	210x180x130mm
Weight	1.5kg

For more information
get in touch with us on +44 (0)1939 234 222