Hanging Scale KERN CH









# The classic among manual hanging scales - perfect for rapid weighing of higher loads

## **Features**

- 11 With the TÜV certification mark, the scales meet the requirements of the standard EN 13155 (Non-fixed load lifting attachments/Breakage resistance) and EN 61010-1 (Electrical safety)
- Ideal for rapid control in goods-in and goods-out
- Also essential in the private sector to determine the weight of fish, game, fruits, bicycle parts, suitcases etc.
- Hold function: For easy reading of the weighing result, the display can be "frozen" in different ways. Either automatically when the weighing value remains unchanged or manually by pressing the Hold key
- Peak load display (peak hold)
- · Hook made of steel

## **Technical data**

- · LCD display, digit height 11 mm
- 2 Tape measure, extractable, length approx.
- Ready for use: Batteries included, 9 V block, operating time up to 20 h
- Weighing units: kg, lb, N
- Overall dimensions W×D×H 90×30×283 mm
- Net weight approx. 0,35 kg
- Permissible ambient temperature 5 °C/35 °C

## Accessories

• 3 Tare pan with pan beam, Max. load 30 kg, details see Accessories, KERN CH-A01N

STANDARD













1	$\mathcal{Z}$	<b>^</b> -				DAkk:
EXT	UNIT	MOVE	BATT	DMS	1 DAY	+3 DAYS

Model	Weighing capacity	Readability	Options  DAkkS Calibr. Certificate
KERN	[Max] kg	[d] g	DARKS CAIDI. CEI IIICATE  DARKS  KERN
CH 15K20	15	20	963-128H
CH 50K50	50	50	963-128H
CH 50K100	50	100	963-128H

# **BALANCES & TEST SERVICE 2024**

**KERN Pictograms** 





## Internal adjusting

Quick setting up of the balance's accuracy with internal adjusting weight (motordriven)



# Adjusting program CAL

For quick setting up of the balance's accuracy. External adjusting weight required



#### **EasyTouch**

Suitable for the connection, data transmission and control through PC or tablet



#### Memory

Balance memory capacity, e.g. for article data, weighing data, tare weights, PLU etc.



## Alibi memory

Secure, electronic archiving of weighing results, complying with the 2014/31/EU standard.



#### **KERN Universal Port** (KUP)

allows the connection of external KUP interface adapters, e.g. RS-232, RS-485, SB, Bluetooth, WIFI, Analogue, Ethernet etc. for the exchange of data and control commands, without installation effort



# RS-232 Data interface

To connect the balance to a printer, PC or network



## **RS-485 Data interface**

To connect the balance to a printer, PC or other peripherals. Suitable for data transfer over large distances. Network in bus topology is possible



# **USB** Data interface

To connect the balance to a printer, PC or other peripherals



#### Bluetooth\* Data interface

To transfer data from the balance to a printer, PC or other peripherals



# WIFI Data interface

To transfer data from the balance to a printer, PC or other peripherals



# **Control outputs**

(optocoupler, digital I/O) To connect relays, signal lamps, valves, etc.



## Analogue interface

to connect a suitable peripheral device for analogue processing of the measurements



#### Interface for second balance

For direct connection of a second balance



#### **Network interface**

For connecting the scale to an Ethernet network



#### **KERN Communication** Protocol (KCP)

It is a standardized interface command set for KERN balances and other instruments, which allows retrieving and controlling all relevant parameters and functions of the device. KERN devices featuring KCP are thus easily integrated with computers, industrial controllers and other digital systems



# GLP/ISO log intern

The balance displays weight, date and time, independent of a printer connection



## **GLP/ISO log Printer**

With weight, date and time. Only with KERN printers.



#### Piece counting

Reference quantities selectable. Display can be switched from piece to weight



#### Recipe level A

The weights of the recipe ingredients can be added together and the total weight of the recipe can be printed out



#### Recipe level B

Internal memory for complete recipés with name and target value of the recipe ingredients. User guidance through display



# Totalising level A

The weights of similar items can be added together and



the total can be printed out Percentage determination



#### Determining the deviation in % from the target value (100 %)

Weighing units Can be switched to e.g. nonmetric units. See



 $\mathcal{Z}$ 

balance model. Please refer to KERN's website for more details



#### Weighing with tolerance range (Checkweighing)

Upper and lower limiting can be programmed individually, e.g. for sorting and dosing. The process is supported by an audible or visual signal, see the relevant model



#### Hold function

(Animal weighing program) When the weighing conditions are unstable, a stable weight is calculated as an average value



# Protection against dust and water splashes IPxx

The type of protection is shown in the pictogram



#### Suspended weighing Load support with hook on the underside of the

balance



# **Battery operation**

Ready for battery operation. The battery type is specified for each device



#### Rechargeable battery pack

Rechargeable set



#### Universal plug-in power supply

with universal input and optional input socket adapters for A) EU, CH, GB B) EU, CH, GB, US C) EU, CH, GB, US, AUS



#### Plug-in power supply 230V/50Hz in standard version for EU, CH.

On request GB, USA or AUS version available



#### Integrated power supply unit

Integrated in balance. 230V/50Hz standard EU. More standards e.g. GB, USA or AUS on request



#### Weighing principle Strain gauges

Electrical resistor on an elastic deforming body



#### Weighing principle Tuning fork

A resonating body is electromagnetically excited, causing it to oscillate



# Weighing principle Electromagnetic force compensation

Coil inside a permanent magnet. For the most accurate weighings



#### Weighing principle Single cell technology

Advanced version of the force compensation principle with the highest level of precision



## Conformity Assessment

The time required for conformity assessment is specified in the pictogram



#### **DAkkS** calibration possible (DKD)

. The time required for DAkkS calibration is shown in days in the pictogram



# Factory calibration (ISO)

The time required for Factory calibration is shown in days in the pictogram



## Package shipment

The time required for internal shipping preparations is shown in days in the pictogram



#### Pallet shipment

The time required for internal shipping preparations is shown in days in the pictogram



<sup>\*</sup>The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by KERN & SOHN GmbH is under license. Other trademarks and trade names are those of their respective owners