Counting Scale KERN CFS







Professional, high-resolution counting scale with 100 item memories and second balance interface, counting resolution up to 75,000 points

Features

- · Memory (PLU) for 100 items with additional text, reference weight and tare weight, e.g. of a container
- Precise counting: The automatic reference weight optimisation of reference weight gradually improves the average piece weight value
- Programmable using numerical key pad:
- required reference quantity
- known reference weight
- · Three displays for weight display, reference weight, total pieces
- Weighing with tolerance range (checkweighing): a visual and audible signal helps with portioning, dispensing or grading
- Fill-to-target function: Programmable target quantity or target weight e.g. for checkweighing. A signal will be displayed when the target value
- · PRE-TARE function for manual subtraction of a known container weight, useful for checking fill-levels

- · Second balance interface to construct a highresolution counting system, standard, e.g. with weighing bridges of KERN KFP V20
- · Draught shield standard for models with weighing plate size A, weighing space W×D×H 155×141×80 mm
- · Protective working cover included with delivery, for models with weighing plate size A, B

Technical data

- · Large backlit LCD displays, digit height 20 mm
- · Dimensions weighing surface
- B W×D 295×225 mm
- C W×D 370×240 mm
- Overall dimensions W×D×H
 - A 315×350×110 mm / 315×350×185 mm (incl. draught shield)
 - **B** 315×350×100 mm
 - C 370×360×100 mm
- Permissible ambient temperature 0 °C/40 °C

Accessories

- · A, B: Protective working cover, scope of delivery: 5 items, KERN CFS-A02S05
- · Internal rechargeable battery pack, operating time up to 70 h without backlight, charging time approx. 14 h, KERN GAB-A04
- Signal lamp for visual support of weighing with tolerance range, KERN CFS-A03
- Y-cable for parallel connection of two terminal devices to the RS-232 interface on the scale, e.g. signal lamp and printer, KERN CFS-A04
- · Further details, plenty of further accessories and suitable printers see Accessories

STANDARI)
	Ī



























OPTION	
	DAkkS
ACCU	+3 DAYS

Model	Weighing capa- city [Max]	Readability [d]	Smallest part weight (Normal)	Counting resolution Points	Net weight approx.	Weighing plate	Options DAkkS Calibr. Certificate DAkkS
KEKIN	kg	g	g/piece	Points	kg		KERN
CFS 300-3	0,3	0,001	0,05	60.000	2,6	A	963-127
CFS 3K-5	3	0,01	0,5	60.000	3,4	В	963-127
CFS 6K0.1	6	0,1	1	60.000	3,2	В	963-128
CFS 15K0.2	15	0,2	2	75.000	3,4	В	963-128
CFS 30K0.5	30	0,5	5	60.000	3,4	В	963-128
CFS 50K-3	50	1	10	50.000	4,4	C	963-128

BALANCES & TEST SERVICE 2024

KERN Pictograms





Internal adjusting

Quick setting up of the balance's accuracy with internal adjusting weight (motordriven)



Adjusting program CAL

For quick setting up of the balance's accuracy. External adjusting weight required



EasyTouch

Suitable for the connection, data transmission and control through PC or tablet



Memory

Balance memory capacity, e.g. for article data, weighing data, tare weights, PLU etc.



Alibi memory

Secure, electronic archiving of weighing results, complying with the 2014/31/EU standard.



KERN Universal Port (KUP)

allows the connection of external KUP interface adapters, e.g. RS-232, RS-485, SB, Bluetooth, WIFI, Analogue, Ethernet etc. for the exchange of data and control commands, without installation effort



RS-232 Data interface

To connect the balance to a printer, PC or network



RS-485 Data interface

To connect the balance to a printer, PC or other peripherals. Suitable for data transfer over large distances. Network in bus topology is possible



USB Data interface

To connect the balance to a printer, PC or other peripherals



Bluetooth* Data interface

To transfer data from the balance to a printer, PC or other peripherals



WIFI Data interface

To transfer data from the balance to a printer, PC or other peripherals



Control outputs

(optocoupler, digital I/O) To connect relays, signal lamps, valves, etc.



Analogue interface

to connect a suitable peripheral device for analogue processing of the measurements



Interface for second balance

For direct connection of a second balance



Network interface

For connecting the scale to an Ethernet network



KERN Communication Protocol (KCP)

It is a standardized interface command set for KERN balances and other instruments, which allows retrieving and controlling all relevant parameters and functions of the device. KERN devices featuring KCP are thus easily integrated with computers, industrial controllers and other digital systems



GLP/ISO log intern

The balance displays weight, date and time, independent of a printer connection



GLP/ISO log Printer

With weight, date and time. Only with KERN printers.



Piece counting

Reference quantities selectable. Display can be switched from piece to weight



Recipe level A

The weights of the recipe ingredients can be added together and the total weight of the recipe can be printed out



Recipe level B

Internal memory for complete recipés with name and target value of the recipe ingredients. User guidance through display



Totalising level A

The weights of similar items can be added together and



the total can be printed out Percentage determination



Determining the deviation in % from the target value (100 %)

Weighing units Can be switched to e.g. nonmetric units. See



 \mathcal{Z}

balance model. Please refer to KERN's website for more details



Weighing with tolerance range (Checkweighing)

Upper and lower limiting can be programmed individually, e.g. for sorting and dosing. The process is supported by an audible or visual signal, see the relevant model



Hold function

(Animal weighing program) When the weighing conditions are unstable, a stable weight is calculated as an average value



Protection against dust and water splashes IPxx

The type of protection is shown in the pictogram



Suspended weighing Load support with hook on the underside of the

balance



Battery operation

Ready for battery operation. The battery type is specified for each device



Rechargeable battery pack

Rechargeable set



Universal plug-in power supply

with universal input and optional input socket adapters for A) EU, CH, GB B) EU, CH, GB, US C) EU, CH, GB, US, AUS



Plug-in power supply 230V/50Hz in standard version for EU, CH.

On request GB, USA or AUS version available



Integrated power supply unit

Integrated in balance. 230V/50Hz standard EU. More standards e.g. GB, USA or AUS on request



Weighing principle Strain gauges

Electrical resistor on an elastic deforming body



Weighing principle Tuning fork

A resonating body is electromagnetically excited, causing it to oscillate



Weighing principle Electromagnetic force compensation

Coil inside a permanent magnet. For the most accurate weighings



Weighing principle Single cell technology

Advanced version of the force compensation principle with the highest level of precision



Conformity Assessment

The time required for conformity assessment is specified in the pictogram



DAkkS calibration possible (DKD)

. The time required for DAkkS calibration is shown in days in the pictogram



Factory calibration (ISO)

The time required for Factory calibration is shown in days in the pictogram



Package shipment

The time required for internal shipping preparations is shown in days in the pictogram



Pallet shipment

The time required for internal shipping preparations is shown in days in the pictogram



^{*}The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by KERN & SOHN GmbH is under license. Other trademarks and trade names are those of their respective owners